

RESPIRATORY
PROTECTION

PAPR Refresher Training

APPLIES TO: All areas where PAPRs are used

How may I obtain refresher training for powered air purifying respirators (PAPRs)?

- Powered Air Purifying Respirator (PAPR) refresher training is available online through the **UC Learning Center**, search for MEDICAL CENTER PAPR TRAINING (1.1) or type in PAPR in the search field.
- Additional resources including PAPR demonstrations over video can be found here:
<https://safety.ucsf.edu/respiratory-protection-program-resources>.
- Each isolation cart stocked with PAPR equipment contains a laminated 'how to use' instructional sheet for the PAPR equipment.
- PAPR training is also available during [fit test](#) sessions for N95 respirators. You do not need to be fit tested to wear a PAPR, however medical clearance to wear any respirator must be on file with Occupational Health Services.
- You may also request instructor-led PAPR training via MCSS Help Desk under 'Safety Training.'

PAPR supplies:

If any supplies are missing for your PAPR or you need to replace the battery, order from Material Services.

If any part of the PAPR is broken, not working or needs repair:

Red tag and immediately remove from service any broken PAPR equipment, file an incident report under Equipment or Environment of Care. Order a replacement PAPR from Material Services.

Fig.	Product	PMM Number (order from Material Services)
A	PAPR Helmet	50332
A	Face shield with chin cuff S/M	50334
A	Face shield with chin cuff M/L	50335
B	Front headband comfort strip	50336
C	Battery (one size)	21920
D	Power cord	235571

Disassembly and Cleaning guidance:

Discard the face shield and comfort strip (if soiled). Safely disconnect and wipe the power cord and battery as well as exposed surfaces of the PAPR helmet between uses.

Face shield securely affixed at each knob across brim



A

Chin cuff up flush in front of ears



B

Power cord should CLICK into place; press tab to release.



C



D